


GNB3 Rabbit mAb

ABZ2658


Product Overview

Name	GNB3 Rabbit mAb
Accession(Primary)	P16520
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (P16520)

Gene Name GNB3**About Gene Symbol**

Western blot analysis of GNB3 in mouse heart
lysates using GNB3 antibody.



Western blot analysis of GNB3 in C6, rat Brain, MCF-7 lysates using GNB3 antibody.