

MGMT Rabbit mAb

ABZ2659


Product Overview

Name	MGMT Rabbit mAb
Accession(Primary)	P16455
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000


Target information (P16455)

Gene Name MGMT

About Gene Symbol



Western blot analysis of MGMT in Jurkat lysates using MGMT antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using MGMT antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.