


MGMT Rabbit mAb

ABZ2660

Product Overview

Name	MGMT Rabbit mAb
Accession(Primary)	P16455
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (P16455)

Gene Name MGMT**About Gene Symbol**

Western blot analysis of MGMT in A549, MCF-7, 293T lysates using MGMT antibody.

