

EpCAM Rabbit mAb


ABZ2662

Product Overview


Name	EpCAM Rabbit mAb
Accession(Primary)	P16422
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35,39 kDa
Recommended Dilution	WB-1000-1:10000; IHC-1:50-1:2000; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (P16422)


Gene Name	EPCAM
About Gene Symbol	




Western blot analysis of EpCAM in rat colon, mouse colon lysates using EpCAM antibody.



Western blot analysis of EpCAM in MCF-7, HCT116, T47D, HeLa lysates using EpCAM antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using EpCAM antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of EpCAM (green) in MCF-7 using EpCAM antibody, and DAPI (blue).