

EpCAM Rabbit mAb


ABZ2663

Product Overview


Name	EpCAM Rabbit mAb
Accession(Primary)	P16422
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 35 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:500-1:5000; FC-1:200-1:500; IP-1:20-1:50

Target information (P16422)


Gene Name	EPCAM
About Gene Symbol	



Western blot analysis of EpCAM in K562 lysates using EpCAM antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using EpCAM antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of EpCAM (green)
in MCF-7 using EpCAM antibody, and DAPI (blue)