

Histone H1.2 Rabbit mAb


ABZ2664

Product Overview


Name	Histone H1.2 Rabbit mAb
Accession(Primary)	P16403
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200

Target information (P16403)

Gene Name	H1-2
About Gene Symbol	



Western blot analysis of Histone H1.2 in mouse lung lysates using Histone H1.2 antibody.



Western blot analysis of Histone H1.2 in HeLa, C6, MCF-7 lysates using Histone H1.2 antibody.