

ALDH2 Rabbit mAb


ABZ2666

Product Overview


Name	ALDH2 Rabbit mAb
Accession(Primary)	P05091
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; IP-1:10-1:100

Target information (P05091)


Gene Name	ALDH2
About Gene Symbol	



Western blot analysis of ALDH2 in mouse liver lysates using ALDH2 antibody.



Western blot analysis of ALDH2 in K562 lysates using ALDH2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using ALDH2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.