

Amyloid Precursor Protein Rabbit mAb


ABZ2668

Product Overview


Name	Amyloid Precursor Protein Rabbit mAb
Accession(Primary)	P05067
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; FC-1:10-1:100; ICC/IF-1:50-1:100

Target information (P05067)


Gene Name	APP
About Gene Symbol	APP (AID(50)) is a cell surface marker involved in cell surface signaling. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of Amyloid Precursor Protein in mouse brain lysates using Amyloid Precursor Protein antibody.



Western blot analysis of Amyloid Precursor Protein in rat Brain lysates using Amyloid Precursor Protein antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Amyloid Precursor Protein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Amyloid Precursor Protein(green) in Hela using Amyloid Precursor Protein antibody, and DAPI(blue).

