

Amyloid Precursor Protein Rabbit mAb


ABZ2669

Product Overview

Name	Amyloid Precursor Protein Rabbit mAb
Accession(Primary)	P05067
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 87 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P05067)

Gene Name	APP
About Gene Symbol	APP (AID(50)) is a cell surface marker involved in cell surface signaling. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Amyloid Precursor Protein in HeLa lysates using Amyloid Precursor Protein antibody.