

Hsp27 Rabbit mAb


ABZ2670

Product Overview


Name	Hsp27 Rabbit mAb
Accession(Primary)	P04792
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:10-1:100; IHC-1:50-1:100

Target information (P04792)


Gene Name	HSPB1
About Gene Symbol	



Western blot analysis of Hsp27 in mouse muscle, mouse heart lysates using Hsp27 antibody.



Western blot analysis of Hsp27 in K562 lysates using Hsp27 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Hsp27 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.