

p53 Rabbit mAb


ABZ2671

Product Overview


Name	p53 Rabbit mAb
Accession(Primary)	P04637
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 44 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500

Target information (P04637)

Gene Name	TP53
About Gene Symbol	TP53 (Tumor Protein P53) is a tumor suppressor protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of p53 in Jurkat lysates using p53 antibody.



Immunocytochemistry analysis of p53 (green) in MCF-7 using p53 antibody, and DAPI (blue)