

p53 Rabbit mAb


ABZ2672

Product Overview

Name	p53 Rabbit mAb
Accession(Primary)	P04637
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 44 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:20-1:50; ICC/IF-1:50-1:100; FC-1:200-1:2000; IP-1:20


Target information (P04637)

Gene Name	TP53
About Gene Symbol	TP53 (Tumor Protein P53) is a tumor suppressor protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of p53 in A431 lysates using p53 antibody.


Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of p53 in T47D lysates using p53 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using p53 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of p53 in HeLa using p53 antibody.