

Von Willebrand Factor Rabbit mAb

ABZ2673


Product Overview

Name	Von Willebrand Factor Rabbit mAb
Accession(Primary)	P04275
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 309 kDa; Observed MW: 309 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:1000-1:5000

Target information (P04275)

Gene Name VWF

About Gene Symbol



← Von Willebrand Factor

Western blot analysis of Von Willebrand Factor in rat kidney, rat lung lysates using von Willebrand Factor antibody.

