


CD90 Rabbit mAb

ABZ2675

Product Overview

Name	CD90 Rabbit mAb
Accession(Primary)	P04216
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P04216)

Gene Name THY1**About Gene Symbol**

Western blot analysis of CD90/Thy1 in rat Brain lysates using CD90 antibody.

