

GBA Rabbit mAb

ABZ2677


Product Overview

Name	GBA Rabbit mAb
Accession(Primary)	P04062
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P04062)


Gene Name GBA

About Gene Symbol



Western blot analysis of GBA in 293T,U87-MG lysates using GBA antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of GBA in U87-MG lysates
using GBA antibody.