

Raf1 Rabbit mAb


ABZ2678

Product Overview


Name	Raf1 Rabbit mAb
Accession(Primary)	P04049
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P04049)

Gene Name	RAF1
About Gene Symbol	RAF1 (Raf-1 Proto-Oncogene, Serine/Threonine Kinase) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Raf1 in C6, 3T3, HeLa lysates using Raf1 antibody.



Immunocytochemistry analysis of Raf1 (green) in MCF-7 using Raf1 antibody, and DAPI (blue)