

HMGCR Rabbit mAb

ABZ2679


Product Overview

Name	HMGCR Rabbit mAb
Accession(Primary)	P04035
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 97 kDa; Observed MW: 97 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (P04035)

Gene Name HMGCR

About Gene Symbol



Western blot analysis of HMGCR in K562, 3T3, HeLa lysates using HMGCR antibody.

