

Vitronectin Rabbit mAb


ABZ2680

Product Overview

Name	Vitronectin Rabbit mAb
Accession(Primary)	P04004
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 54 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100

Target information (P04004)

Gene Name	VTN
About Gene Symbol	



Western blot analysis of Vitronectin in Human serum lysates using Vitronectin antibody.