

MMP1 Rabbit mAb

ABZ2681


Product Overview

Name	MMP1 Rabbit mAb
Accession(Primary)	P03956
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1; 50-1:100; IP-1:20

Target information (P03956)

Gene Name MMP1

About Gene Symbol



Copyright© 2024 Zen-bioscience Limited.

