

Transferrin Rabbit mAb


ABZ2684

Product Overview

Name	Transferrin Rabbit mAb
Accession(Primary)	P02787
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; FC-1:20-1:50

Target information (P02787)

Gene Name	TF
About Gene Symbol	



Western blot analysis of Transferrin in K562
lysates using Transferrin antibody.