

Transferrin Receptor Rabbit mAb


ABZ2686

Product Overview


Name	Transferrin Receptor Rabbit mAb
Accession(Primary)	P02786
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (P02786)

Gene Name	TFRC
About Gene Symbol	TFRC (Transferrin Receptor) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Transferrin Receptor in K562, HeLa lysates using Transferrin Receptor 1 antibody.



Immunocytochemistry analysis of Transferrin Receptor(green) in HeLa using Transferrin Receptor antibody, and DAPI(blue)