


PF4 Rabbit mAb

ABZ2687


Product Overview

Name	PF4 Rabbit mAb
Accession(Primary)	P02776
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200


Target information (P02776)

Gene Name PF4**About Gene Symbol**

Western blot analysis of PF4 in rat spleen lysates using PF4 antibody.



Western blot analysis of PF4 in A549, HL-60 lysates using PF4 antibody.



Immunohistochemistry analysis of paraffin-embedded liver cancer using PF4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.