

Prealbumin Rabbit mAb

ABZ2688


Product Overview

Name	Prealbumin Rabbit mAb
Accession(Primary)	P02766
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:500

Target information (P02766)

Gene Name TTR

About Gene Symbol



Western blot analysis of Prealbumin in rat Brain lysates using Prealbumin antibody.

