

MBP Rabbit mAb

ABZ2689


Product Overview

Name	MBP Rabbit mAb
Accession(Primary)	P02686
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 33 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:200-1:500

Target information (P02686)

Gene Name MBP

About Gene Symbol



Western blot analysis of Myelin Basic Protein in rat Brain lysates using Myelin Basic Protein antibody.

