

MBP Rabbit mAb

ABZ2690


Product Overview

Name	MBP Rabbit mAb
Accession(Primary)	P02686
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P02686)

Gene Name MBP

About Gene Symbol



Western blot analysis of Myelin Basic Protein in HeLa, CHO-K1, C6 lysates using Myelin Basic Protein antibody.

