



Collagen I alpha 1 Rabbit mAb

ABZ2694


Product Overview

Name	Collagen I alpha 1 Rabbit mAb
Accession(Primary)	P02452
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 139 kDa; Observed MW: 139 kDa
Recommended Dilution	WB-1:1000-1:2000; IF: 1:50; IP-1:20-1:50

Target information (P02452)

Gene Name COL1A1**About Gene Symbol**

Immunocytochemistry analysis of Collagen I alpha 1 (green) in Hela using Collagen I alpha 1 antibody, and DAPI(blue).



Western blot analysis of Collagen I in A549, HL-60, U2OS lysates using Collagen I antibody.