

CD8 alpha Rabbit mAb


ABZ2695

Product Overview


Name	CD8 alpha Rabbit mAb
Accession(Primary)	P01732
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 26 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; IHC-1:50-1:100

Target information (P01732)


Gene Name	CD8A
About Gene Symbol	



Western blot analysis of CD8 in THP-1 lysates using CD8 alpha antibody.



Western blot analysis of CD8 in K562, CHO-K1, HeLa lysates using CD8 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CD8 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.