

## TGF alpha Rabbit mAb


ABZ2696

### Product Overview


|                             |                                                                         |
|-----------------------------|-------------------------------------------------------------------------|
| <b>Name</b>                 | TGF alpha Rabbit mAb                                                    |
| <b>Accession(Primary)</b>   | P01135                                                                  |
| <b>Host Species</b>         | Rabbit                                                                  |
| <b>Clonality</b>            | Monoclonal                                                              |
| <b>Isotype</b>              | IgG                                                                     |
| <b>Species Reactivity</b>   | Human                                                                   |
| <b>Calculated MW</b>        | Calculated MW: 17 kDa; Observed MW: 17 kDa                              |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:500-1:2000;<br>ICC/IF-1:50-1:100; IP-1:20-1:50 |

### Target information (P01135)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | TGFA |
| <b>About Gene Symbol</b> |      |



Western blot analysis of TGF alpha in K562 lysates using TGF alpha antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using TGF alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.