


PDGF B Rabbit mAb

ABZ2697

Product Overview

Name	PDGF B Rabbit mAb
Accession(Primary)	P01127
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:50-1:100

Target information (P01127)

Gene Name PDGFB**About Gene Symbol**

Western blot analysis of PDGF B in mouse brain, rat brain lysates using PDGF B antibody.

