

GTPase NRAS Rabbit mAb


ABZ2698

Product Overview

Name	GTPase NRAS Rabbit mAb
Accession(Primary)	P01111
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:10-1:100

Target information (P01111)

Gene Name	NRAS
About Gene Symbol	



Western blot analysis of Ras in HeLa, 3T3, C6, rat Brain, MCF-7 lysates using GTPase NRAS antibody.