


GTPase HRAS Rabbit mAb

ABZ2699


Product Overview

Name	GTPase HRAS Rabbit mAb
Accession(Primary)	P01111
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P01111)

Gene Name NRAS**About Gene Symbol**

Western blot analysis of GTPase HRAS in mouse brain lysates using GTPase HRAS antibody.



Western blot analysis of RAS in rat Brain lysates using RAS antibody.