

GTPase HRAS Rabbit mAb


ABZ2700

Product Overview


Name	GTPase HRAS Rabbit mAb
Accession(Primary)	P01111
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:200-1:500; FC-1:20; IP-1:20

Target information (P01111)

Gene Name	NRAS
About Gene Symbol	



Western blot analysis of GTPase HRAS in 3T3, mouse brain lysates using GTPase HRAS antibody.



Western blot analysis of Ras in HeLa, CHO-K1, C6 lysates using Ras antibody.