

c-Myc Rabbit mAb


ABZ2701

Product Overview


Name	c-Myc Rabbit mAb
Accession(Primary)	P01106
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 57-65 kDa
Recommended Dilution	WB-1:1000-1:4000; IHC-1:200-1:800; IF: 1:50; FC-1:50-1:100; IP-1:20-1:50

Target information (P01106)

Gene Name	MYC
About Gene Symbol	MYC (MYC Proto-Oncogene, BHLH Transcription Factor) is a proto-oncogene involved in translational regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of c-Myc (green) in HeLa using c-Myc antibody, and DAPI (blue).



Western blot analysis of c-Myc in K562, C6, 3T3, HeLa lysates using c-Myc antibody.