

Cystatin C Rabbit mAb


ABZ2702

Product Overview


Name	Cystatin C Rabbit mAb
Accession(Primary)	P01034
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 16 kDa; Observed MW: 13 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P01034)


Gene Name	CST3
About Gene Symbol	CST3 (Cystatin C) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cystatin C in rat lung lysates using Cystatin C antibody.



Western blot analysis of Cystatin C in K562, 3T3 lysates using Cystatin C antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Cystatin C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.