

C3 Rabbit mAb


ABZ2703

Product Overview

Name	C3 Rabbit mAb
Accession(Primary)	P01024
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 187 kDa; Observed MW: 187 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-:1:200-1:2000

Target information (P01024)

Gene Name	C3
About Gene Symbol	C3 (Complement C3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of C3 in Human serum
lysates using C3 antibody.