

## C3 Rabbit mAb


ABZ2704

### Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | C3 Rabbit mAb                                |
| <b>Accession(Primary)</b>   | P01024                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human                                        |
| <b>Calculated MW</b>        | Calculated MW: 187 kDa; Observed MW: 120 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; ICC/IF-1:50-1:100; FC-1:20  |

### Target information (P01024)

|                          |                                                                                                                                                                                                                                                                                                     |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | C3                                                                                                                                                                                                                                                                                                  |
| <b>About Gene Symbol</b> | C3 (Complement C3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of C3 in K562, HeLa lysates using C3 antibody.