


alpha 1 Antitrypsin Rabbit mAb

ABZ2705

Product Overview

Name	alpha 1 Antitrypsin Rabbit mAb
Accession(Primary)	P01009
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 47 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (P01009)

Gene Name SERPINA1**About Gene Symbol**

Western blot analysis of alpha 1 Antitrypsin in HeLa, A549, MCF-7, U251, HL-60 lysates using alpha 1-Antitrypsin antibody.

