

ASS1 Rabbit mAb


ABZ2707

Product Overview


Name	ASS1 Rabbit mAb
Accession(Primary)	P00966
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:10-1:100

Target information (P00966)


Gene Name	ASS1
About Gene Symbol	




Western blot analysis of ASS1 in rat liver lysates using ASS1 antibody.



Western blot analysis of ASS1 in HeLa lysates using ASS1 antibody.



Immunocytochemistry analysis of ASS1(green) in Hela using ASS1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using ASS1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.