

PGK1 Rabbit mAb


ABZ2708

Product Overview


Name	PGK1 Rabbit mAb
Accession(Primary)	P00558
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P00558)


Gene Name	PGK1
About Gene Symbol	PGK1 (Phosphoglycerate Kinase 1) is a protein kinase and signaling regulator involved in cell migration and invasion. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of PGK1 in C6, Raw264.7, PC-12, 3T3 lysates using PGK1 antibody.



Western blot analysis of PGK1 in HeLa lysates using PGK1 antibody.



Western blot analysis of PGK1 in MCF-7, HCT116, T47D, HeLa lysates using PGK1 antibody



Immunocytochemistry analysis of PGK1(green) in HeLa using PGK1 antibody, and DAPI(blue).