

Aspartate Aminotransferase Rabbit mAb


ABZ2709

Product Overview


Name	Aspartate Aminotransferase Rabbit mAb
Accession(Primary)	P00505
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; IP-1:10-1:100

Target information (P00505)

Gene Name	GOT2
About Gene Symbol	GOT2 (Glutamic-Oxaloacetic Transaminase 2) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and



Western blot analysis of FABP1 in K562, rat Brain, 3T3 lysates using Aspartate Aminotransferase antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using FABP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.