

Superoxide Dismutase 1 Rabbit mAb

ABZ2711


Product Overview

Name	Superoxide Dismutase 1 Rabbit mAb
Accession(Primary)	P00441
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:20; ICC/IF-1:20; FC-1:20


Target information (P00441)

Gene Name SOD1

About Gene Symbol



Western blot analysis of Superoxide Dismutase 1 in K562, rat Brain, C6, 3T3, HeLa lysates using Superoxide Dismutase 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Superoxide Dismutase 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.