

MTCO2 Rabbit mAb


ABZ2714

Product Overview


Name	MTCO2 Rabbit mAb
Accession(Primary)	P00403
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 26 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (P00403)


Gene Name	MT-CO2
About Gene Symbol	



Western blot analysis of MTCO2 in K562, HeLa lysates using MTCO2 antibody.



Immunocytochemistry analysis of MTCO2(green) in HeLa using MTCO2 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MTCO2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.