


MTCO2 Rabbit mAb

ABZ2715


Product Overview

Name	MTCO2 Rabbit mAb
Accession(Primary)	P00403
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 26 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (P00403)

Gene Name MT-CO2**About Gene Symbol**

Western blot analysis of MTCO2 in K562, C6, HeLa lysates using MTCO2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MTCO2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.