

ALDH1A1 Rabbit mAb


ABZ2718

Product Overview


Name	ALDH1A1 Rabbit mAb
Accession(Primary)	P00352
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20-1:50

Target information (P00352)

Gene Name	ALDH1A1
About Gene Symbol	



Western blot analysis of ALDH1A1 in K562, 3T3 lysates using ALDH1A1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using ALDH1A1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.