

Lactate Dehydrogenase A Rabbit mAb


ABZ2719

Product Overview

Name	Lactate Dehydrogenase A Rabbit mAb
Accession(Primary)	P00338
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:20-1:100

Target information (P00338)

Gene Name	LDHA
About Gene Symbol	



Western blot analysis of LDHA in Jurkat, rat Brain, C6, HeLa lysates using Lactate Dehydrogenase A antibody.