

Lactate Dehydrogenase A Rabbit mAb

ABZ2720


Product Overview


Name	Lactate Dehydrogenase A Rabbit mAb
Accession(Primary)	P00338
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:5000-1:150000; IHC-1:200-1:2000; IF: 1:50

Target information (P00338)


Gene Name LDHA

About Gene Symbol





Immunocytochemistry analysis of Lactate Dehydrogenase(green) in HeLa using Lactate Dehydrogenase antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Lactate Dehydrogenase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.