

Cyclin E2 Rabbit mAb


ABZ2721

Product Overview


Name	Cyclin E2 Rabbit mAb
Accession(Primary)	O96020
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:100; IHC-1:50-1:100

Target information (O96020)


Gene Name	CCNE2
About Gene Symbol	CCNE2 (Cyclin E2) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Cyclin E2 in Jurkat lysates using Cyclin E2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Cyclin E2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Cyclin E2
(green) in MCF-7 using Cyclin E2 antibody, and
DAPI(blue)