

BAG3 Rabbit mAb


ABZ2727

Product Overview

Name	BAG3 Rabbit mAb
Accession(Primary)	O95817
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 62 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:50-1:100

Target information (O95817)

Gene Name	BAG3
About Gene Symbol	BAG3 (BAG Cochaperone 3) is a molecular chaperone involved in apoptosis. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Bag3 in C6, 3T3, HeLa lysates using BAG3 antibody.