


Epac1 Rabbit mAb


ABZ2728

Product Overview

Name	Epac1 Rabbit mAb
Accession(Primary)	O95398
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 104 kDa; Observed MW: 104 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (O95398)

Gene Name RAPGEF3**About Gene Symbol**



Western blot analysis of Epac1 in 3T3 lysates
using Epac1 antibody.