


INMT Rabbit mAb

ABZ2731


Product Overview

Name	INMT Rabbit mAb
Accession(Primary)	O95050
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:100-1:200

Target information (O95050)

Gene Name INMT**About Gene Symbol**

Western blot analysis of INMT in SKBR3, mouse lung lysates using INMT antibody.



Western blot analysis of INMT in C6 lysates using
INMT antibody.